


XOP – STATUS



ROGER DEJUS
Control Account Manager for
Insertion Devices for APS Upgrade

Outline

- History and Scope
- Current distribution (v2.4)
- Overview
- Examples
- References
- Summary

History and Scope


- Developed during the early 1990s to suit local needs at the ESRF and the APS (efforts officially merged 1995; now 25 years in the making ... end of life)
- XOP v2.0: 1,000 CD-ROMs distributed (400 registered users)
- Front-end graphical user interface for computer codes of different origins and different languages for the synchrotron radiation community (written in the licensed Interactive Data Language IDL)
 - Modelling of x-ray sources
 - Characterization of optical elements (mirrors, filters, crystals, multilayers, etc.) and their effect on sources (“pipes” via files)
 - Multipurpose visualizations and data analyses
 - Optional plug-in of external software packages “extensions” expands the functionality of XOP

Current XOP Distribution v2.4 (2015)



- Used worldwide by synchrotron radiation facilities (and many others) and has been crucial for beamline designers and users for over a decade. The current version runs on Unix, Linux, Mac OS X, and Windows.
- The graphical user interface and many modules of the code are written in IDL, which is subject to U.S. Export Control and is categorized under Export Control Classification Number (ECCN) 5D002 (as of IDL v8.0). As such, it can only be distributed to users who have completed and submitted an application that is approved by the Argonne's Export Control process.
- Licensed since April 2015 with 2,029 application requests to date (10% denied)
- Licensed with IDL v8.3 embedded and good thru September 20, 2021 (no additional IDL embedding requests planned)
- Replaced by the **OASYS** software, which contains similar functionalities (codes) but also with enhancements and new codes for different applications

Overview

X-Ray Sources




X-Ray Optics and Photon-Atom Interactions



XOP Extensions


- SHADOWVUI
 - Interface to the SHADOW ray-tracing code (Cerrina et al.)
- IMD
 - Multilayer software (Windt)
- TOPO
 - Surface topography (Windt)

General Purpose Tools and Documentation





XPLOT


XOP Examples – XTC and XINPRO




On-Axis
Brilliance Tuning
Curves (XTC)





Si X-Ray Crystal
Diffraction Profile
(XINPRO)



XOP Examples – XUS: Helical SCU vs. APS Undulator A




Power Density h-SCU left,
UA right



UA Horizontal
Linear
Polarization

- $K_y = 1.635$
- Distance 26.5 m
- General appearance correct – sharp edges not




h-SCU Circular
Polarization

- $K_x = K_y = 1.20$
- Distance 26.5 m
- General appearance correct – sharp edges not

7

3rd Harmonic (18 keV); very low flux (not important)

XOP Request for Downloads



Software Download for XOP v2.4 and Future Releases

Software Description

XOP (X-ray Oriented Programs) is a widget-based computer program used as a common front-end interface of codes of interest to the synchrotron radiation community. It provides codes formodeling of x-ray sources (e.g., synchrotron radiation sources, such as undulators and wigglers), characteristics of optical elements(mirror, filters, crystals, multilayers, etc.), and multipurpose data visualizations and analyses.

Instructions

New users please select the "Apply" button and fill out the application form.
Authorized users please select the "Login" button and use your credentials to login.
After Export Control Review you will be notified by e-mail how to use the "Authorized User Login."

Terms of Use

The software's intended use is to calculate radiation properties of x-ray sources and their propagation through optical elements.This software is subject to U.S. Export Control because it uses the Interactive Data Language (IDL) from Exelis Visual, and it is categorized under Export Control Classification Number (ECCN) 5D002. The software applicant agrees that the software and documentation will not be shipped, transferred or exported into any country, or used in any manner prohibited by U.S. export laws as detailed herein.

A successful registration and Export Control screening allows the software applicant to operate, install (including on-site installation for multiple user access), maintain, repair, overhaul and refurbishing version 2.4 of this software and future releases of the same unless the ECCN is changed.

Licensee acknowledges that it is not a citizen, national, or resident of, and is not under the control of the government of Cuba, Iran, North Korea, Sudan or Syria, and that Licensee will not download or otherwise export or re-export the Software Product and any related technical data directly or indirectly to the above mentioned countries nor to citizens, nationals, or residents of those countries. Licensee will comply with all applicable export, re-export and foreign policy controls and restrictions imposed by the United States and will take the necessary actions and precautions to ensure that it complies with all such laws or regulations.

- I am not located in and that I am not a national of Cuba, Iran, North Korea, Syria, or Sudan, and I agree not to transfer the Software Product to any such locations or individuals;
- I am not listed or affiliated with anyone on any of the following lists:
 - [Denied Persons List](#) - U.S. Department of Commerce
 - [Unverified List](#) - U.S. Department of Commerce
 - [Entity List](#) - Supplement No. 4 to part 744 of the Export Administration Regulations
 - [Specially Designated Nationals](#) - U.S. Treasury
- I also agree to comply with all U.S. and international laws (including the obtaining of any required import license or authorization), and I understand and agree that Exelis Visual Information Solutions is not responsible for my failure to comply with any such laws.

If not self-employed, you hereby certify that you are authorized to act on behalf of your company.

View Exelis Visual Information Solutions (Harris Geospatial Solutions) [Export Compliance](#) document for information regarding the export of the Software Product. We recommend you seek legal counsel if you have any questions regarding the legality of your use of the Software Product.

The applicant needs to accept these terms on the XOP Software Download Application Form.

Start a New Application

Apply

Authorized User Login

Login

- <https://beam.aps.anl.gov/apps/xop/>

References

- Manuel Sánchez del Río and Roger J. Dejus "XOP v2.4: recent developments of the x-ray optics software toolkit", Proc. SPIE 8141, Advances in Computational Methods for X-Ray Optics II, 814115 (23 September 2011); <https://doi.org/10.1117/12.893911>
- <https://www.aps.anl.gov/Science/Scientific-Software>
- <https://beam.aps.anl.gov/apps/xop/>
- srio@esrf.fr or srio@lbl.gov
- dejus@anl.gov

XOP SUCCESSFUL AND APPRECIATED 25 YEARS DOWNLOAD REQUESTS CONTINUES REPLACED BY THE NEW OASYS SOFTWARE



**U.S. DEPARTMENT OF
ENERGY**

Argonne National Laboratory is a
U.S. Department of Energy laboratory
managed by UChicago Argonne, LLC.



An aerial photograph of the Argonne National Laboratory campus, showing various buildings, parking lots, and surrounding greenery, all overlaid with a semi-transparent blue filter.

THANK YOU!

QUESTIONS?



U.S. DEPARTMENT OF
ENERGY

Argonne National Laboratory is a
U.S. Department of Energy laboratory
managed by UChicago Argonne, LLC.

Argonne 
NATIONAL LABORATORY